



Dynamic Diagnostics, Inc.

Material Safety Data Sheet

Section 1. Chemical Product and Company Information

Common Name:	95% Ethyl Alcohol (Denatured)	Code:	9521
Supplier:	Dynamic Diagnostics, Inc. 12665 Richfield Court Livonia, MI 48150 800-717-4677	MSDS#:	9521
Synonym:	None.	Validation Date:	09-01-2008
Trade Name:	None.	Print Date:	09-01-2008
Material Uses:	Dehydrating solution.	Responsible Name:	Dr. B
Manufacturer:	BBC Biochemical 409 Eleanor Lane Mount Vernon, WA 98273	In Case of Emergency: 1-800-424-9300 Chemtrec USA 1-202-483-7616 Chemtrec Intrl 1-800-635-4477	

Section 2. Composition and Information on Ingredients

Name	CAS#	OSHA PEL/TWA	Exposure Limits
1) Denatured Ethyl Alcohol	64-17-5		None established
2) Water	7732-18-5		

Section 3. Hazards Identification

Physical State and Appearance	Colorless liquid with characteristic alcoholic odor.
Emergency Overview	Not available.
Routes of Entry	Inhalation, ingestion, and skin absorption.
Potential Acute Health Effects	
Eyes	Irritation.
Skin	Irritation.
Inhalation	May cause dizziness, headache, nausea, and narcosis.
Ingestion	May cause blindness, nausea, damage to gastrointestinal tract, liver, kidneys, and cardiovascular system.
Potential Chronic Health Effects	Components not listed (IARC, NTP, OSHA) as cancer causing agents.
Medical Conditions Aggravated by Overexposure/Overexposure/Signs/Symptoms	Kidney, liver, heart, and GI conditions.

Section 4. First Aid Measures

Eye Contact	Immediately flush thoroughly with water for at least 15 minutes.
Skin Contact	Wash thoroughly with soap and water.
Inhalation	Remove to fresh air; give artificial respiration if breathing has stopped.
Ingestion	If conscious, drink water and induce vomiting.
Notes to Physician	Not available.

Section 5. Fire Fighting Measures

Flammability of the	Not available.
---------------------	----------------

Product	
Auto-ignition Temperature	Not available.
Flash Points	About 55° F
Flammable Limits	Not available.
Products of Combustion	Not available.
Fire Hazards in Presence Of Various Substances	Not available.
Explosion Hazards in Presence of Various Substances	Not available.
Fire Fighting Media and Instructions	Dry chemical, carbon dioxide, alcohol foam. Use water spray to cool fire-exposed containers and disperse vapors.
Protective Clothing (Fire)	Fire fighters should use self-contained breathing apparatus and protective clothing.
Special Remarks on Fire Hazards	Vapor can travel distances to ignition source and flash back.
Special Remarks on Explosion Hazards	

Section 6. Accidental Release Measures

Small Spill and Leak	Evacuate the area of all unnecessary personnel. Wear suitable protective equipment. Eliminate all sources of ignition and provide ventilation.
Large Spill and Leak	Same as above.

Section 7, Handling and Storage

Handling	Use adequate ventilation. Avoid repeated or prolonged skin contact or breathing of vapors. Do not breathe vapors. Do not get in eyes.
Storage	Store and use in areas away from heat, sparks, and open flame. Store in tightly closed containers in a cool, dry, well-ventilated, fire-resistant area. Store away from oxidizing agents.

Section 8. Exposure Controls / Personal Protection

Engineering Controls	
Personal Protection	
Eyes	Safety glasses with approved side shields should be worn at all times.
Body	Protective clothing, such as a lab coat, to minimize skin contact.
Respiratory	Use NIOSH or MSHA approved respirators in absence of proper environmental control.
Hands	Neoprene, Nitrile or equivalent protective gloves.
Feet	Impervious footwear should be worn to prevent skin contact.
Personal Protection in Case of a Large Spill	See above.
Product Name	Exposure Limits
1)	
2)	
3)	
Consult Local authorities before acceptable exposure limits.	

Section 9. Physical and Chemical Properties

Physical State and Appearance	Colorless liquid.	Odor:	Characteristic alcoholic odor.
Molecular Weight	Not applicable.	Taste:	Vinous.
Molecular Formula	Not applicable.	Color:	Colorless.
pH (1%/Water)	Not applicable.		
Boiling/Condensation Point	173° F		
Melting/Freezing Point	Not applicable.		
Critical Temperature	Not applicable.		
Specific Gravity	Not applicable.		
Vapor Pressure	Not applicable.		
Vapor Density	Not applicable.		
Volatility	Not applicable.		
Odor Threshold	Not applicable.		
Evaporation Rate	Not applicable.		
VOC	Not available.		
Viscosity	Not available.		
Ionicity (in water)	Not available.		

Dispersion Properties	Not available.
Solubility	Complete in water.
Physical Chemical Comments	Not available.

Section 10. Stability and Reactivity

Stability and Reactivity	Stable.
Conditions of Instability	Heat.
Incompatibility with Various Substances	Oxidizers.
Hazardous Decomposition Products	Thermal-oxidative degradation can produce oxides of carbon. Toxic gases and vapors (i.e. carbon monoxide, formaldehyde) may be released in a fire.
Hazardous Polymerization	Does not occur.

Section 11. Toxicological Information

Toxicity to Animals	LD50: Not available. LC50: Not available.
Chronic Effects on Humans	This product contains the following toxic chemical subject to the reporting requirements of section 313 of the Emergency Planning and Community Right-To-Know Act of 1986 and of 40 CFR 372: CAS # _____ Chemical Name _____ DeMinimis/SARA 313% NONE
Other Toxic Effects on Humans	
Special Remarks on Toxicity to Animals	Not available.
Special Remarks on Chronic Effects on Humans	Not available.
Special Remarks on Other Toxic Effects on Humans	Not available.

Section 12. Ecological Information

Ecotoxicity	Not available.
BODS and COD	Not available.
Biodegradable/OEDC	Not available.
Mobility	Not available.
Toxicity of the Products of Biodegradation	Not available.
Special Remarks on The Products of Biodegradation	Not available.

Section 13. Disposal Considerations

Waste Information	Not available.
Waste Stream	Not available.
Consult your local or regional authorities.	

Section 14. Transport Information

DOT Classification	Flammable liquid.
Marine Pollutant	Not available.
Hazardous Substances Reportable Quantity	Not available.
Special Provisions for Transport	Not applicable.
TDG Classification	Not controlled under TDG (Canada).
ADR/RID Classification	Not controlled under ADR (Europe).
IMO/IMDG Classification	Not controlled under IMDG.
ICAO/IATA Classification	Not controlled under IATA.


Section 15. Other Information

Label requirements

Hazardous Material Information System (U.S.A.)

Health	1
Fire Hazard	3
Reactivity	0
Personal Protection	

National Fire Protection Association (U.S.A.)



References

Other Special Considerations

Notice to Reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.